

Lipofectamine 3000 protocol

PLASMID / CONDITION	NG/ML	WELLS	
EV	823	4.5	×
WT	718	4.5	×
9 wells total			

1. Seed cells in a 12 well plate

- Cell type: Caco-2

2. No. of cells: _____ (typical $0.5\text{--}2 \times 10^5$)

3. Cells should be 70–90 % confluent.

4. Work in the dark

5. Reduce the medium volume in each well to **0.5 mL** to improve transfection efficiency.

6. Prepare two mixes

a. Optimem and Lipofectamine (official range 1.5–3 μL per μg DNA, this sheet uses the low end)

- 50 μL Optimem per well 450 μL
- 1.5 $\mu\text{L}/\mu\text{g}$ DNA $\rightarrow$ **0.75 μL** Lipofectamine per well 6.75 μL

b. Optimem, DNA and P3000 Reagent – add P3000 last (change the DNA and Lipofectamine and P3000 follow)

- 50 μL Optimem per well 450 μL
- 500 ng DNA per well 4500 ng
- 2 $\mu\text{L}/\mu\text{g}$ DNA $\rightarrow$ **1 μL** P3000 per well 9 μL

7. Mix well (no need to vortex)

8. Add DNA Mix to the Lipofectamine mix

9. Incubate 10–15 min and drop-wise **100 μL** mix per well.

10. Check cells in microscope

11. Spin down plate

12. ~~Change medium at the end of the day (2.5 h minimum)~~

13. Incubate cells for 2–4 days at 37 °C

EV 4.5 wells

 225 μL Optimem
3.375 μL Lipofectamine

 225 μL Optimem
2.73 μL DNA
4.5 μL P3000

WT 4.5 wells

 225 μL Optimem
3.375 μL Lipofectamine

 225 μL Optimem
3.13 μL DNA
4.5 μL P3000

EV: $500 \text{ ng} \times 4.5 \div 823 \text{ ng}/\mu\text{L} = 2.73 \mu\text{L}$

WT: $500 \text{ ng} \times 4.5 \div 718 \text{ ng}/\mu\text{L} = 3.13 \mu\text{L}$