

cDNA synthesis

29 Sept 2026

1000 ng RNA per sample · RNA + water to 15.00 µL · then 4.00 µL reaction mix + 1.00 µL reverse transcriptase · 20.00 µL per tube

#	sample	ng/µL	water µL	RNA µL
1	Ev 1	300.559	11.67 ☒	3.33 ☒
2	Ev 2	449.094	12.77 ☒	2.23 ☐
3	Wt 1	373.196	12.32 ☒	2.68 ☐
4	Wt 2	372.41	12.31 ☒	2.69 ☐
5	Mut1	536.236	13.14 ☐	1.86 ☐
6	Mut 2	687.963	13.55 ☐	1.45 ☐
7	Ev ap 1	249.556	10.99 ☐	4.01 ☐
8	Ev ap 2	415.826	12.60 ☐	2.40 ☐
9	Wt ap1	554.331	13.20 ☐	1.80 ☐
10	Wt ap 2	805.519	13.76 ☐	1.24 ☐
11	Mut ap 1	485.669	12.94 ☐	2.06 ☐
12	Mut ap 2	341.141	12.07 ☐	2.93 ☐

Master mix × 14 (12 samples + 10 %)

Reaction mix 4.00 × 14 56.00 µL ☐

Reverse transcriptase 1.00 × 14 14.00 µL ☐

Add 5.00 µL to each tube.

