

cDNA synthesis

08 Oct 2026

1000 ng RNA per sample · RNA + water to 15.00 µL · then 4.00 µL reaction mix + 1.00 µL reverse transcriptase · 20.00 µL per tube

# sample	ng/µL	water µL	RNA µL
1 Ev 1	300.56	11.67 ☒	3.33 ☐
2 Ev 2	449.09	12.77 ☐	2.23 ☐
3 Wt 1	373.20	12.32 ☐	2.68 ☐
4 Wt 2	372.41	12.31 ☐	2.69 ☐
5 Mut1	536.24	13.14 ☐	1.86 ☐
6 Mut 2	687.96	13.55 ☐	1.45 ☐
7 Ev ap 1	249.56	10.99 ☐	4.01 ☐
8 Ev ap 2	415.83	12.60 ☐	2.40 ☐
9 Wt ap1	554.33	13.20 ☐	1.80 ☐
10 Wt ap 2	805.52	13.76 ☐	1.24 ☐
11 Mut ap 1	485.67	12.94 ☐	2.06 ☐
12 Mut ap 2	341.14	12.07 ☐	2.93 ☐

Master mix × 14 (12 samples + 10 %)

Reaction mix 4.00 × 14 56.00 µL ☐

Reverse transcriptase 1.00 × 14 14.00 µL ☐

Add 5.00 µL to each tube.

